



HEINE BETA® 200 Ophthalmoscope

The ophthalmoscope with superior aspherical optics



The ophthalmoscope with unique aspherical optics and all the features of the BETA200S, except for the macula aperture/pinhole and the following **range of lenses: 27 lenses from -35 D to +40 D.**

+ in 1 D steps:	1 – 10	15	20	40	D
- in 1 D steps:	1 – 10	15	20	25	35 D

BETA200 Ophthalmoscope	2.5V	3.5V
without handle	C-001.30.100	C-002.30.100
BETA200 M2 Ophthalmoscope, without handle. Like BETA200, except with fixation star aperture without polar coordinates and with pinhole instead of small spot.		
	C-001.30.102	C-002.30.102
XHL Xenon Halogen spare bulb	X-001.88.069	X-002.88.070

- ∴ Eliminates corneal reflex
- ∴ Ideal for small pupils
- ∴ Dustproof

BETA200 apertures



Slit, fixation star with polar coordinates, cobalt blue filter, large spot, small spot, hemispot



With red-free filter

BETA200 M2 apertures



Slit, fixation star, cobalt blue filter, large spot, pinhole, hemispot



With red-free filter